

**High-speed steel, 4% chromium, 5% cobalt und 5% molybdenum, no. 1.3243 - chips**

Art. ID                      ECRM-D 290-1  
Unit                         100 g  
Deliverydetails            No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	0,911	%		
	Silicon (Si)	[7440-21-3]	0,072	%		
	Manganese (Mn)	[7439-96-5]	0,244	%		
	Phosphorus (P)	[7723-14-0]	0,016	%		
	Sulfur (S)	[7704-34-9]	0,016	%		
	Chromium (Cr)	[7440-47-3]	4,18	%		
	Molybdenum (Mo)	[7439-98-7]	4,83	%		
	Nickel (Ni)	[7440-02-0]	0,329	%		
	Cobalt (Co)	[7440-48-4]	5,12	%		
	Copper (Cu)	[7440-50-8]	0,081	%		
	Vanadium (V)	[7440-62-2]	1,91	%		
	Tungsten (W)	[7440-33-7]	6,27	%		
	Alloy / Werkstoff		1.3243			