

2,2',3,5'-Tetrabromodiphenyl ether (BDE 44) 10 µg/mL in 2,2,4-Trimethylpentane (Isooctane)

Art. ID LM24-OP-4000869
Unit 10 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2

Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,2',3,5'-Tetrabromodiphenyl ether (BDE 44)	[446254-20-2]	10	µg/mL		