

Aluminium (Al) 1000 mg/L in diluted HCl, Single-Element Standard for ICP

Art. ID LM24-OP-4001023
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Aluminium (Al) | [7429-90-5] | 1000 | mg/L | | |