

Beryllium (Be) 1000 mg/L in 2-5% HCl, Single-Element Standard for ICP

Art. ID LM24-OP-4001053
Unit 30 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1760 III E1

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Beryllium (Be) | [7440-41-7] | 1000 | mg/L | | |