

**Iridium (Ir) 1000 mg/L in 5-10% HCl, Single-Element Standard for ICP**

Art. ID LM24-OP-4001157  
Unit 30 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1789 II E2

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Iridium (Ir)      | [7439-88-5] | 1000                | mg/L |        |        |