

Tin (Sn) 1000 mg/L in 2-5% HNO₃/HF tr, Single-Element Standard for ICP

Art. ID LM24-OP-4001311
Unit 250 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Tin (Sn) | [7440-31-5] | 1000 | mg/L | | |