

**Tin (Sn) 1000 mg/L in 20% HCl, Single-Element Standard for ICP**

Art. ID LM24-OP-4001316

Unit 500 mL

Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1789 II E2

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Tin (Sn)          | [7440-31-5] | 1000                | mg/L |        |        |