

Tellurium (Te) 1000 mg/L in 10-20% HCl, Single-Element Standard for ICP

Art. ID LM24-OP-4001335

Unit 250 mL

Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1789 II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|--------------|---------------------|------|--------|--------|
| | Tellurium (Te) | [13494-80-9] | 1000 | mg/L | | |