

Zirconium (Zr) 1000 mg/L in 2-5% HNO₃/HF tr, Single-Element Standard for ICP

Art. ID LM24-OP-4001401
Unit 30 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN3264 III E1

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Zirconium (Zr) | [7440-67-7] | 1000 | mg/L | | |