

Arsenic (As) 10000 mg/L in 2-5% HNO₃, Single-Element Standard for ICP

Art. ID LM24-OP-4001424
Unit 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Arsenic (As) | [7440-38-2] | 10000 | mg/L | | |