

Copper (Cu) 10000 mg/L in 5% HNO₃, Single-Element Standard for ICP

Art. ID LM24-OP-4001486

Unit 100 mL

Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Copper (Cu) | [7440-50-8] | 10000 | mg/L | | |