

Iron (Fe) - 10 g/L in diluted HCl for ICP

Art. ID LM24-OP-4001504
Unit 250 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG)

Description

Validity: 12 months

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Iron (Fe) | [7439-89-6] | 10 | g/L | | |