

Potassium (K) 10000 mg/L in diluted HNO₃, Single-Element Standard for ICP

Art. ID LM24-OP-4001543
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Potassium (K) | [7440-09-7] | 10000 | mg/L | | |