

Molybdenum (Mo) 10000 mg/L in 2-5 % NH₄OH, Single-Element Standard for ICP

Art. ID LM24-OP-4001577
Unit 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Molybdenum (Mo) | [7439-98-7] | 10000 | mg/L | | |