

**Molybdenum (Mo) 10000 mg/L in 2-5% HNO<sub>3</sub>/HF tr, Single-Element Standard for ICP**

Art. ID LM24-OP-4001578  
Unit 30 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Molybdenum (Mo)	[7439-98-7]	10000	mg/L		