

**Molybdenum (Mo) 10000 mg/L in diluted HNO<sub>3</sub>/HF, Single-Element Standard for ICP**

Art. ID LM24-OP-4001580

Unit 250 mL

Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2922 (6.1) II E2

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Molybdenum (Mo)   | [7439-98-7] | 10000               | mg/L |        |        |