

Tin (Sn) 10000 mg/L in diluted HNO₃/HF tr, Single-Element Standard for ICP

Art. ID LM24-OP-4001690
Unit 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2922 (6.1) II E2

Description

Conservation 36 months (in originally closed bottle).

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Tin (Sn) | [7440-31-5] | 10000 | mg/L | | |