

Strontium (Sr) 10000 mg/L in 2-5% HCl, Single-Element Standard for ICP

Art. ID LM24-OP-4001694
Unit 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Strontium (Sr) | [7440-24-6] | 10000 | mg/L | | |