

Tantalum (Ta) - 10 g/L in diluted HNO3/HF for ICP

Art. ID LM24-OP-4001700
Unit 30 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN2922 (6.1) II E2

Description

Validity: 12 months

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Tantalum (Ta) | [7440-25-7] | 10 | g/L | | |