

Tantalum (Ta) 10000 mg/L in diluted HNO₃/HF Single-Element Standard for ICP

Art. ID LM24-OP-4001701
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2922 (6.1) II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Tantalum (Ta) | [7440-25-7] | 10000 | mg/L | | |