

**Tantalum (Ta) 10000 mg/L in diluted HNO<sub>3</sub>/HF tr, Single-Element Standard for ICP**

Art. ID LM24-OP-4001703  
Unit 500 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2922 (6.1) II E2

**Description**

Conservation 36 months (in originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tantalum (Ta)	[7440-25-7]	10000	mg/L		