

**Titanium (Ti) 10000 mg/L in 2-5% HCl/HF tr, Single-Element Standard for ICP**

Art. ID LM24-OP-4001712  
Unit 30 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1760 III E1

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Titanium (Ti)     | [7440-32-6] | 10000               | mg/L |        |        |