

Tungsten (W) 10000 mg/L in 5-10% NH4OH, Single-Element Standard for ICP

Art. ID LM24-OP-4001740
Unit 30 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ)

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Tungsten (W) | [7440-33-7] | 10000 | mg/L | | |