

Tungsten (W) 10000 mg/L in diluted HNO₃/HF tr, Single-Element Standard for ICP

Art. ID LM24-OP-4001747
Unit 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2922 (6.1) II E2

Description

Conservation 36 months (in originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tungsten (W)	[7440-33-7]	10000	mg/L		