

Zinc (Zn) 10000 mg/L in 5% HNO₃, Single-Element Standard for ICP

Art. ID LM24-OP-4001761
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Zinc (Zn) | [7440-66-6] | 10000 | mg/L | | |