

ICP-MS MIX 10 components in HCl 2% , HF tr

Art. ID LM24-OP-4001777
Unit 100 mL
Deliverydetails No Dangerous Good /not restricted

Description

Shelf-life: 12 months

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|--------------|---------------------|------|--------|--------|
| | Tungsten (W) | [7440-33-7] | 40 | mg/L | | |
| | Rhodium (Rh) | [7440-16-6] | 40 | mg/L | | |
| | Platinum (Pt) | [7440-06-4] | 40 | mg/L | | |
| | Palladium (Pd) | [7440-05-3] | 40 | mg/L | | |
| | Tellurium (Te) | [13494-80-9] | 40 | mg/L | | |
| | Zirconium (Zr) | [7440-67-7] | 40 | mg/L | | |
| | Silicon (Si) | [7440-21-3] | 40 | mg/L | | |
| | Molybdenum (Mo) | [7439-98-7] | 50 | mg/L | | |
| | Antimony (Sb) | [7440-36-0] | 80 | mg/L | | |
| | Tin (Sn) | [7440-31-5] | 80 | mg/L | | |