

ICP-MS MIX 10 components in HCl 2% , HF tr

Art. ID LM24-OP-4001777
Unit 100 mL
Deliverydetails No Dangerous Good /not restricted

Description

Shelf-life: 12 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tungsten (W)	[7440-33-7]	40	mg/L		
	Rhodium (Rh)	[7440-16-6]	40	mg/L		
	Platinum (Pt)	[7440-06-4]	40	mg/L		
	Palladium (Pd)	[7440-05-3]	40	mg/L		
	Tellurium (Te)	[13494-80-9]	40	mg/L		
	Zirconium (Zr)	[7440-67-7]	40	mg/L		
	Silicon (Si)	[7440-21-3]	40	mg/L		
	Molybdenum (Mo)	[7439-98-7]	50	mg/L		
	Antimony (Sb)	[7440-36-0]	80	mg/L		
	Tin (Sn)	[7440-31-5]	80	mg/L		