

Standard Solution 5 components in HNO3 5%

Art. ID LM24-OP-4001782
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Validity: 12 months, Standard equivalent to Jobin Yvon

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Potassium (K)	[7440-09-7]	1500	µg/mL		
	Lead (Pb)	[7439-92-1]	1000	µg/mL		
	Aluminium (Al)	[7429-90-5]	500	µg/mL		
	Magnesium (Mg)	[7439-95-4]	500	µg/mL		
	Cadmium (Cd)	[7440-43-9]	100	µg/mL		