

ICP-MS MIX 4 components in 2% HNO₃, 0.1% HF

Art. ID LM24-OP-4001794
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Description

Shelf-life: 12 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Molybdenum (Mo)	[7439-98-7]	10	mg/L		
	Antimony (Sb)	[7440-36-0]	10	mg/L		
	Tin (Sn)	[7440-31-5]	10	mg/L		
	Titanium (Ti)	[7440-32-6]	10	mg/L		