

**Internal standard mix for ICP-MS systems - 8 components in HNO₃ 10% , HCl tr%
(Standard equivalent to Agilent/ Varian Equivalent)**

Art. ID LM24-OP-4001819
Unit 100 mL
Delivery details ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1789 II E2

Description

Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Lithium-6 (Li-6)	[14258-72-1]	10	µg/mL		
	Scandium (Sc)	[7440-20-2]	10	µg/mL		
	Germanium (Ge)	[7440-56-4]	10	µg/mL		
	Rhodium (Rh)	[7440-16-6]	10	µg/mL		
	Indium (In)	[7440-74-6]	10	µg/mL		
	Terbium (Tb)	[7440-27-9]	10	µg/mL		
	Lutetium (Lu)	[7439-94-3]	10	µg/mL		
	Bismuth (Bi)	[7440-69-9]	10	µg/mL		