

ICP-MS MIX 12 components in 2% HNO₃

Art. ID LM24-OP-4001822
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Shelf-life: 8 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	2,5	mg/L		
	Cadmium (Cd)	[7440-43-9]	0,5	mg/L		
	Chromium (Cr)	[7440-47-3]	5	mg/L		
	Mercury (Hg)	[7439-97-6]	0,2	mg/L		
	Nickel (Ni)	[7440-02-0]	5	mg/L		
	Lead (Pb)	[7439-92-1]	5	mg/L		
	Zinc (Zn)	[7440-66-6]	5	mg/L		
	Barium (Ba)	[7440-39-3]	2,5	mg/L		
	Copper (Cu)	[7440-50-8]	5	mg/L		
	Molybdenum (Mo)	[7439-98-7]	3	mg/L		
	Antimony (Sb)	[7440-36-0]	0,5	mg/L		
	Selenium (Se)	[7782-49-2]	0,7	mg/L		