

Quality Control Standard of Mineral Oils (ISO 9377-2-Mod), 2 components, 500 mg/L of each component in Acetone

Art. ID LM24-OP-4001852
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1090 II E2

Description

Validity: 24 months

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-----------------------------|-------------|---------------------|------|--------|--------|
| | Mineral Oil [CAS:8042-47-5] | [8042-47-5] | 500 | mg/L | | |
| | Diesel Oil | [n/a] | 500 | mg/L | | |