

**Florisil Cartridge QC Standard Mixture (ISO 9377-2-Mod)2 components, 1000 mg/L of each component in Hexane**

Art. ID LM24-OP-4001853  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1208 II E2

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Mineral Oil [CAS:8042-4 7-5]	[8042-47-5]	1000	mg/L		
	Diesel Oil	[n/a]	1000	mg/L		