

Quality Control Standard of Mineral Oils (ISO 9377-2-Mod)2 components, 10000 mg/L of each component in Hexane

Art. ID LM24-OP-4001855
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1208 II E2

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Mineral Oil [CAS:8042-4 7-5]	[8042-47-5]	10000	mg/L		
	Diesel Oil	[n/a]	10000	mg/L		