

**Stearyl Stearate Test Solution (ISO 9377-2-Mod)1 el., Stearyl stearate 2000 µg/mL,  
n-Hexane 1part**

Art. ID LM24-OP-4001859  
Unit 10 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1208 II E2

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Stearyl stearate	[2778-96-3]	2000	µg/mL		