

**n-Alkane Standard Solution - 16 components C10-C40 (all even) (ISO 9377-2-Mod), 100 mg/L of each component in Hexane/Petroleum ether (1/1)**

Art. ID LM24-OP-4001860  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1993 II E2

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	n-Decane (C10)	[124-18-5]	100	mg/L		
	n-Dodecane (C12)	[112-40-3]	100	mg/L		
	n-Tetradecane (C14)	[629-59-4]	100	mg/L		
	n-Hexadecane (C16)	[544-76-3]	100	mg/L		
	n-Octadecane (C18)	[593-45-3]	100	mg/L		
	n-Eicosane (C20)	[112-95-8]	100	mg/L		
	n-Docosane (C22)	[629-97-0]	100	mg/L		
	n-Tetracosane (C24)	[646-31-1]	100	mg/L		
	n-Hexacosane (C26)	[630-01-3]	100	mg/L		
	n-Octacosane (C28)	[630-02-4]	100	mg/L		
	n-Triacontane (C30)	[638-68-6]	100	mg/L		
	n-Dotriacontane (C32)	[544-85-4]	100	mg/L		
	n-Tetratriacontane (34)	[14167-59-0]	100	mg/L		
	n-Hexatriacontane (C36)	[630-06-8]	100	mg/L		
	n-Octatriacontane (C38)	[7194-85-6]	100	mg/L		
	n-Tetracontane (C40)	[4181-95-7]	100	mg/L		