

n-Alkane Standard Solution 16 components C10÷C40 (all even) (ISO 9377-2-Mod) in n-Hexane

Art. ID LM24-OP-4001861
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1208 II E2

Description
Shelf-life: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Decane, No. of C: C10:0	[124-18-5]	50	mg/L		
	n-Dodecane	[112-40-3]	50	mg/L		
	n-Tetradecane	[629-59-4]	50	mg/L		
	Hexadecane, No. of C: C16:0	[544-76-3]	50	mg/L		
	n-Octadecane	[593-45-3]	50	mg/L		
	Eicosane, No. of C: C20:0	[112-95-8]	50	mg/L		
	Docosane, No. of C: C22:0	[629-97-0]	50	mg/L		
	n-Tetracosane	[646-31-1]	50	mg/L		
	Hexacosane, No. of C: C26:0	[630-01-3]	50	mg/L		
	n-Octacosane	[630-02-4]	50	mg/L		
	n-Triacontane	[638-68-6]	50	mg/L		
	Dotriacontane, No. of C: C32:0	[544-85-4]	50	mg/L		
	n-Tetratriacontane	[14167-59-0]	50	mg/L		
	Hexatriacontane, No. of C: C36:0	[630-06-8]	50	mg/L		
	n-Octatriacontane	[7194-85-6]	50	mg/L		
	n-Tetracontane	[4181-95-7]	50	mg/L		