

n-Alkane Standard Solution - 6 components (ISO 9377-2-Mod) 100 µg/mL of each component in Hexane

Art. ID LM24-OP-4001863
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1208 II E2

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	n-Decane (C10)	[124-18-5]	100	µg/mL		
	n-Hexadecane (C16)	[544-76-3]	100	µg/mL		
	n-Eicosane (C20)	[112-95-8]	100	µg/mL		
	n-Tetracosane (C24)	[646-31-1]	100	µg/mL		
	n-Triacontane (C30)	[638-68-6]	100	µg/mL		
	n-Tetracontane (C40)	[4181-95-7]	100	µg/mL		