

BTEX Standard Solution - 6 components (ISO 9377-2-Mod), 200 mg/L of each component in Hexane

Art. ID LM24-OP-4001864
Unit 1 mL
Deliverydetails De Minimis UN1208 II E2

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Benzene	[71-43-2]	200	mg/L		
	Ethylbenzene	[100-41-4]	200	mg/L		
	Toluene	[108-88-3]	200	mg/L		
	o-Xylene	[95-47-6]	200	mg/L		
	m-Xylene	[108-38-3]	200	mg/L		
	p-Xylene	[106-42-3]	200	mg/L		