

Internal Stock Standard Solution - for aromatic hydrocarbons 2 components (ISO 15009), 2000 µg/mL of each component in Methanol

Art. ID LM24-OP-4001875
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Toluene D8	[2037-26-5]	2000	µg/mL		
	Ethylbenzene-D10	[25837-05-2]	2000	µg/mL		