

**Stock calibration standard Solution - 4 - 22 components (ISO 15680,2003), 2000 µg/mL of each component in Methanol**

Art. ID LM24-OP-4001883  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description  
Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Cumene	[98-82-8]	2000	µg/mL		
	Bromobenzene	[108-86-1]	2000	µg/mL		
	1,1,2,2-Tetrabromoethane	[79-27-6]	2000	µg/mL		
	1,2,3-Trichloropropane	[96-18-4]	2000	µg/mL		
	n-Propylbenzene	[103-65-1]	2000	µg/mL		
	2-Chlorotoluene	[95-49-8]	2000	µg/mL		
	4-Chlorotoluene	[106-43-4]	2000	µg/mL		
	1,3,5-Trimethylbenzene	[108-67-8]	2000	µg/mL		
	tert.-Butylbenzene	[98-06-6]	2000	µg/mL		
	1,2,4-Trimethylbenzene	[95-63-6]	2000	µg/mL		
	sec-Butylbenzene	[135-98-8]	2000	µg/mL		
	1,3-Dichlorobenzene	[541-73-1]	2000	µg/mL		
	1,4-Dichlorobenzene	[106-46-7]	2000	µg/mL		
	4-Isopropyltoluene	[99-87-6]	2000	µg/mL		
	1,2-Dichlorobenzene	[95-50-1]	2000	µg/mL		
	n-Butylbenzene	[104-51-8]	2000	µg/mL		
	1,2-Dibromo-3-chloropropane	[96-12-8]	2000	µg/mL		
	1,2,4-Trichlorobenzene	[120-82-1]	2000	µg/mL		
	Hexachloro-1,3-butadiene	[87-68-3]	2000	µg/mL		
	Naphthalene	[91-20-3]	2000	µg/mL		
	1,2,3-Trichlorobenzene	[87-61-6]	2000	µg/mL		
	1,3,5-Trichlorobenzene	[108-70-3]	2000	µg/mL		