

## Daily Quality Control Standard (without Internal Standard) (ASTM D5580-02) - 14 components

Art. ID LM24-OP-4001912  
Unit 1 mL  
Deliverydetails De Minimis

### Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	n-Hexane (C6)	[110-54-3]	12	%		
	n-Heptane (C7)	[142-82-5]	20	%		
	n-Octane (C8)	[111-65-9]	15	%		
	n-Decane (C10)	[124-18-5]	10	%		
	n-Dodecane (C12)	[112-40-3]	1	%		
	2,2,4-Trimethylpentane (Isooctane)	[540-84-1]	20	%		
	Benzene	[71-43-2]	1	%		
	Toluene	[108-88-3]	9	%		
	Ethylbenzene	[100-41-4]	2	%		
	p-Xylene	[106-42-3]	3	%		
	o-Xylene	[95-47-6]	2	%		
	1,2,4-Trimethylbenzene	[95-63-6]	3	%		
	1,2,4,5-Tetramethylbenz ene	[95-93-2]	1	%		
	Naphthalene	[91-20-3]	1	%		