

Daily Quality Control Standard (with Internal Standard) (ASTM D5580-02) - 15 components

Art. ID LM24-OP-4001913
Unit 1 mL
Deliverydetails De Minimis

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	n-Hexane (C6)	[110-54-3]	11	%		
	n-Heptane (C7)	[142-82-5]	19	%		
	n-Octane (C8)	[111-65-9]	14	%		
	n-Decane (C10)	[124-18-5]	9	%		
	n-Tridecane (C13)	[629-50-5]	1	%		
	2,2,4-Trimethylpentane (Isooctane)	[540-84-1]	19	%		
	Benzene	[71-43-2]	1	%		
	Toluene	[108-88-3]	8	%		
	Ethylbenzene	[100-41-4]	1	%		
	p-Xylene	[106-42-3]	2	%		
	o-Xylene	[95-47-6]	1	%		
	1,2,4-Trimethylbenzene	[95-63-6]	2	%		
	1,2,4,5-Tetramethylbenz ene	[95-93-2]	1	%		
	Naphthalene	[91-20-3]	1	%		
	2-Hexanone	[591-78-6]	10	%		