

ACs Standard Solution - 25 components 50 mg/L of each component in n-Pentane

Art. ID LM24-OP-4001926
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1265 II E2

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Benzene	[71-43-2]	50	mg/L		
	Toluene	[108-88-3]	50	mg/L		
	Ethylbenzene	[100-41-4]	50	mg/L		
	o-Xylene	[95-47-6]	50	mg/L		
	m-Xylene	[108-38-3]	50	mg/L		
	p-Xylene	[106-42-3]	50	mg/L		
	n-Heptane (C7)	[142-82-5]	50	mg/L		
	n-Octane (C8)	[111-65-9]	50	mg/L		
	n-Nonane (C9)	[111-84-2]	50	mg/L		
	n-Decane (C10)	[124-18-5]	50	mg/L		
	n-Dodecane (C12)	[112-40-3]	50	mg/L		
	n-Tetradecane (C14)	[629-59-4]	50	mg/L		
	n-Hexadecane (C16)	[544-76-3]	50	mg/L		
	n-Octadecane (C18)	[593-45-3]	50	mg/L		
	n-Eicosane (C20)	[112-95-8]	50	mg/L		
	n-Docosane (C22)	[629-97-0]	50	mg/L		
	n-Tetracosane (C24)	[646-31-1]	50	mg/L		
	n-Hexacosane (C26)	[630-01-3]	50	mg/L		
	n-Octacosane (C28)	[630-02-4]	50	mg/L		
	n-Triacontane (C30)	[638-68-6]	50	mg/L		
	n-Dotriacontane (C32)	[544-85-4]	50	mg/L		
	n-Tetracontane (C34)	[14167-59-0]	50	mg/L		
	n-Hexatriacontane (C36)	[630-06-8]	50	mg/L		
	n-Octatriacontane (C38)	[7194-85-6]	50	mg/L		
	n-Tetracontane (C40)	[4181-95-7]	50	mg/L		