

PAH Standard Solution - 7 components, 10 µg/mL of each component in Toluene

Art. ID LM24-OP-4001939
Unit 1 mL
Deliverydetails De Minimis UN1294 II E2

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Benzo(a)anthracene	[56-55-3]	10	µg/mL		
	Benzo(b)fluoranthene	[205-99-2]	10	µg/mL		
	Benzo(j)fluoranthene	[205-82-3]	10	µg/mL		
	Benzo(k)fluoranthene	[207-08-9]	10	µg/mL		
	Indeno(1,2,3-cd)pyrene	[193-39-5]	10	µg/mL		
	Dibenz(a,h)anthracene	[53-70-3]	10	µg/mL		
	Benzo(a)pyrene	[50-32-8]	10	µg/mL		