

PAH Standard Solution - 4 components, 100 µg/mL of each component in Acetonitrile

Art. ID LM24-OP-4001940
Unit 1 mL
Deliverydetails De Minimis UN1648 II E2

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Benzo(g,h,i)perylene D1 2	[93951-66-7]	100	µg/mL		
	Benzo(a)pyrene D12	[63466-71-7]	100	µg/mL		
	Naphthalene D8	[1146-65-2]	100	µg/mL		
	Fluoranthene D10	[93951-69-0]	100	µg/mL		