

## **PCBs Standard Solution - 8 components, 1 µg/mL of each component in Hexane**

Art. ID LM24-OP-4001944  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1208 II E2

### Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,2',3,3',4,4',5,5'-Octachlorobiphenyl (PCB 194)	[35694-08-7]	1	µg/mL		
	2,2',3,4,4',5,5'-Heptachlorobiphenyl (PCB 180)	[35065-29-3]	1	µg/mL		
	2,2',4,4',5,5'-Hexachlorobiphenyl (PCB 153)	[35065-27-1]	1	µg/mL		
	2,2',3,4,4',5'-Hexachlorobiphenyl (PCB 138)	[35065-28-2]	1	µg/mL		
	2,3',4,4',5'-Pentachlorobiphenyl (PCB 118)	[31508-00-6]	1	µg/mL		
	2,2',4,5,5'-Pentachlorobiphenyl (PCB 101)	[37680-73-2]	1	µg/mL		
	2,2',5,5'-Tetrachlorobiphenyl (PCB 52)	[35693-99-3]	1	µg/mL		
	2,4,4'-Trichlorobiphenyl (PCB 28)	[7012-37-5]	1	µg/mL		