

OPP Internal Standard Solution - 1 component (EPA 507) Triphenyl phosphate 500 µg/mL in Methyl-tert.butylether

Art. ID LM24-OP-4001988
Unit 1 mL
Deliverydetails De Minimis UN2398 II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|---------------------|------------|---------------------|-------|--------|--------|
| | Triphenyl phosphate | [115-86-6] | 500 | µg/mL | | |