

Decachlorobiphenyl Standard Solution 1 component (EPA 508.A)
2,2',3,3',4,4',5,5',6,6'-Decachlorobiphenyl (PCB 209) 1000 µg/mL in Toluene

Art. ID LM24-OP-4002017
Unit 1 mL
Deliverydetails De Minimis UN1294 II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------------------------------------------|-------------|---------------------|-------|--------|--------|
| | 2,2',3,3',4,4',5,5',6,6'-Decachlorobiphenyl (PCB 209) | [2051-24-3] | 1000 | µg/mL | | |