

**PAH Standard Solution - 13 components (EPA 525.1), 500 µg/mL of each component in Toluene**

Art. ID LM24-OP-4002041  
Unit 1 mL  
Deliverydetails De Minimis UN1294 II E2

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Acenaphthylene	[208-96-8]	500	µg/mL		
	Anthracene	[120-12-7]	500	µg/mL		
	Benzo(a)anthracene	[56-55-3]	500	µg/mL		
	Benzo(j)fluoranthene	[205-82-3]	500	µg/mL		
	Benzo(k)fluoranthene	[207-08-9]	500	µg/mL		
	Benzo(g,h,i)perylene	[191-24-2]	500	µg/mL		
	Benzo(a)pyrene	[50-32-8]	500	µg/mL		
	Chrysene	[218-01-9]	500	µg/mL		
	Dibenz(a,h)anthracene	[53-70-3]	500	µg/mL		
	Fluorene	[86-73-7]	500	µg/mL		
	Indeno(1,2,3-cd)pyrene	[193-39-5]	500	µg/mL		
	Phenanthrene	[85-01-8]	500	µg/mL		
	Pyrene	[129-00-0]	500	µg/mL		